



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

TMEM125-transmembrane protein 125 Gene

Reference: PTXBC016858

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	6330530A05Rik; AV223516; transmembrane protein 125		
For research use only.			

GENE DESCRIPTION	transmembrane protein 125
GENE ID	128218
ORIGIN SPECIES	Human
ACCESSIONS	BC016858
VECTOR	pDONR223
NCBI SYMBOL	TMEM125
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]