



Reference: PTXBC016860

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAIG-3; RAIG3; G-protein coupled receptor family C group 5 member C; orphan G-protein coupled receptor; retinoic acid responsive gene protein; retinoic acid-induced gene 3 protein; G protein-coupled receptor class C group 5 member C		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor, family C, group 5, member C
GENE ID	55890
ORIGIN SPECIES	Human
ACCESSIONS	BC016860
VECTOR	pDONR223
NCBI SYMBOL	GPRC5C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr