

## ZCCHC4-zinc finger, CCHC domain containing 4 Gene

Reference: PTXBC016914

|                |                                                                                                                                                                                                                         |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | HSPC052; ZGRF4; zinc finger CCHC domain-containing protein 4; DHHC domain-containing zinc finger protein; zinc finger, CCHC domain containing 4; zinc finger, GRF-type containing 4; zinc finger CCHC-type containing 4 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | zinc finger, CCHC domain containing 4                                                             |
| <b>GENE ID</b>          | 29063                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC016914                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ZCCHC4                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggaggtggtctctcttggatcctgccgtccccgcccgctgtgccctcaaggaccactctctgtttgtaaaggtgacccaagggaagaagaactcggaggtttatgcctgttcagcctgtagatagaaaaagactgtaatttttcagtggaagatgaaagttgtcaggagct  
agacttgcgtcccgagaagctcataaccgaagatgtcagctccccgtgtccgaacgcagtggtgtggaaggtacttgaagttattgagttgcccttgactcagagaaagtgttgcacacatgtcagcagttgtgtaccagatgactggggccaacatagtgagcatcagttctgggtaa  
tgtgtccattaccagttgaagaaggcccgactcaactctttatccactggaaaacaagaagaatgccagtatctgtttgctgatcggagctgtcagttctgttagacttacttctgccctcggaatcagaagagtactgtgtgtggaaccaaggttgcagtgatgaagttgaca  
gcatcaggtgacaagaagtctaattaaaagcctttattggatattgatttgcgtag