



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ANGEL2-angel homolog 2 (Drosophila) Gene

Reference: PTXBC016966

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Ccr4d; KIAA0759L; protein angel homolog 2; angel homolog 2		
For research use only.			

GENE DESCRIPTION	angel homolog 2 (Drosophila)
GENE ID	90806
ORIGIN SPECIES	Human
ACCESSIONS	BC016966
VECTOR	pDONR223
NCBI SYMBOL	ANGEL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]