



Reference: PTXBC016967

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf19; C18B11; RNA pseudouridylylase domain-containing protein 2; C18B11 homolog (44.9kD); RNA pseudouridylylase domain containing 2		

For research use only.

GENE DESCRIPTION	RNA pseudouridylate synthase domain containing 2
GENE ID	27079
ORIGIN SPECIES	Human
ACCESSIONS	BC016967
VECTOR	pDONR223
NCBI SYMBOL	RPUSD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr