



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

EPHA4-EPH receptor A4 Gene

Reference: PTXBC016981

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEK8; SEK; TYRO1; ephrin type-A receptor 4; EK8; EPH-like kinase 8; TYRO1 protein tyrosine kinase; receptor protein-tyrosine kinase HEK8; tyrosine-protein kinase TYRO1; tyrosine-protein kinase receptor SEK; EPH receptor A4		

For research use only.

GENE DESCRIPTION	EPH receptor A4
GENE ID	2043
ORIGIN SPECIES	Human
ACCESSIONS	BC016981
VECTOR	pDONR223
NCBI SYMBOL	EPHA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]