



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NECAP2-NECAP endocytosis associated 2 Gene

Reference: PTXBC017014

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	adaptin ear-binding coat-associated protein 2; adaptin-ear-binding coat-associated protein 2; NECAP endocytosis associated 2		

For research use only.

GENE DESCRIPTION	NECAP endocytosis associated 2
GENE ID	55707
ORIGIN SPECIES	Human
ACCESSIONS	BC017014
VECTOR	pDONR223
NCBI SYMBOL	NECAP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagagagagcgggtacgagtcggtgctgctggtgctcaagcgtacgctacacgtctaccgcatccctccggcggtaccacacgtggctacagggctcgaggtgagcgtgaccacccatcagagtgccgctgaggtactcgtcaaaaggacagatggcctacatcaagctg
gaggacagagcgtcaggggagcctcttctcaggccccggtgatcaattctctggcacagcgtgtggagaggtgtagcggattccagcaggtactctgtgatccgcctcgaagatggaattggcgacggcgctttattggaattggctctggggaccgaggtgatgccttgcactcaatgt
tgcattgcaggacatttcaagtgggtgaacaacagatgtgaatttgcaaaacaaagccgaacccagacacaaagccctaaatcgaccctgggctctcaagaggggccagacatcaagctcgaacatcgaacaatgaagaagaagaaggaagcagctgggaatcccgagtcggcc
tgccagcacagaggggctgagcctgcttccccctccccaggggggaaaacctccacctgatccctccccctggggagcaggtgtgtagcgtgggatccctctccagacacagtggtgccagcttgcaggtctgga
ctggaggacatttaccaaatcagagatcaacttcacagcagacccagcagcggctggctgcaggtctga