

DNAJC12-DnaJ (Hsp40) homolog, subfamily C, member 12 Gene

Reference: PTXBC017018

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dnaJ homolog subfamily C member 12; DnaJ (Hsp40) homolog, subfamily C, member 12; J domain containing protein 1 (JDP1); J domain protein 1; j domain-containing protein 1; DnaJ heat shock protein family (Hsp40) member C12		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 12
GENE ID	56521
ORIGIN SPECIES	Human
ACCESSIONS	BC017018
VECTOR	pDONR223
NCBI SYMBOL	DNAJC12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgcaatactgaattacaggtcagaagatactgaagattactacacattactgggatgtgatgaactatcttgggtgaacaaatcctggcagaatttaaagtcagagctctggaatgtcaccagacaagcatcctgaaaaccccaagctgtggagactttcagaactgcagaagg
caaggagattctgaccaatgaagagagtcgagcccgctatgaccactggcgaaggagccagatgtcgtccattccagcagtgaggaaagctttgaatgactcagtgaaagacgtcaatgcactgggtgtcagaggtaaaaaagacgtgatctggaagaatctgacaagactcatacca
ccaagatggaaaatgaggaatgtaatgagcaaaagagaagaaagaaagaggagctggcttcaaccgagagaaaacggagcagaagaacccaagccctagagaagtcagctctcccgcataatcagattctcaggttttcagatgtgaatgggtggcacctctgttcgcgtggt
ccaaggatgctcctcagaactcctgaggaagtcagaaactatgaaatata