

UNC50-unc-50 homolog (C. elegans) Gene

Reference: PTXBC017019

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GMH1; HSD23; PDLs22; URP; protein unc-50 homolog; geal-6 membrane-associated high-copy suppressor 1; periodontal ligament-specific protein 22; protein GMH1 homolog; unc-50 related; unc-50-like protein; uncoordinated-like protein; unc-50 inner nuclear memb		

For research use only.

GENE DESCRIPTION	unc-50 homolog (C. elegans)
GENE ID	25972
ORIGIN SPECIES	Human
ACCESSIONS	BC017019
VECTOR	pDONR223
NCBI SYMBOL	UNC50
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtaccgagtactcagtgaaattccttagtcaggggaacggagcttgaattccagggaatgcggcaagacacacagccggagcgaacgctacaataatctgagaaggctttccgcttccgcaaatggacttgaatttgcctctggcagatgctctacctgtcacatccccacaga
gagtttacagaaatttcattatcgaacacagacgaaggaccagtgggccagagatgacctgcttcttggctctgttaagtatctggctctgtgtccactataggatttggcttctgtctggacatgggattctttgagacaataaagcttctcttgggtgtactcatagattgttaggcgtt
ggctctctgatagcaactttaatgtgttcatctctaacaagtatttagtgaacgacagacagagactatgatgtggaatggggctatgcttttgatgtgcatctcaatgcttttatccactcctggctatttgcattttatccagctttttcatcaaccatgttatcctgacagacacattatggat
atttagttggaataaccttatgggttggtgcagttggctattatatactatgaacttctctgggatacagtgcatggcatttttgaataacagtaattctctgtatcatttgcacctctgattctgctctacgggctttccctggcactgggatggaactcaccatactctgttcttctataaglac
agagtgaataa