

GJB2-gap junction protein, beta 2, 26kDa Gene

Reference: PTXBC017048

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CX26; DFNA3; DFNA3A; DFNBI; DFNBI1; HID; KID; NSRD1; PPK; gap junction beta-2 protein; connexin 26; gap junction protein, beta 2, 26kDa; mutant gap junction protein beta 2; gap junction protein beta 2		

For research use only.

GENE DESCRIPTION	gap junction protein, beta 2, 26kDa
GENE ID	2706
ORIGIN SPECIES	Human
ACCESSIONS	BC017048
VECTOR	pDONR223
NCBI SYMBOL	GJB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattggggcagcgtgcagacgatcctgggggtgtgaacaacactccaccagcattgaaagatctggtcaccgtctcttctatttgcattatgatcctcgttggctgcaaggaggtgtggggagatgagcaggccgactttgtctgcaacacctgcagccaggctgcaa
gaacgtgtgtacgatcactacttcccacatcccacatccggctatgggcccctgcagctgatcttctgttccacgccagcgtcctcctagtgccatgcacgtggcctaccggagacatgagaagaagaggaagtcatcaaggggagataaagagtgaatttaaggacatcaggagat
caaaaccagaaggtccgatcgaaggtccctgtgggtgacacacaaagcagcatcttctccgggtcatctcgaagccgcctcatgacgtcttctatgcatgacgacgcttctccatgcagcggctgtgtgaagtgcacgcctggcctgtcccaactgtggactgttgtgt
cccgcccccaggagaagactgttccagtggtcatgattgcagtgctgtgaattgcatcctgctgaatgtcactgaattgtgtattgtctaattagatattgttgggaagtcacaaaagccagtttaa