

CHCHD4-coiled-coil-helix-coiled-coil-helix domain containing 4 Gene

Reference: PTXBC017082

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TIMM40; mitochondrial intermembrane space import and assembly protein 40; coiled-coil-helix-coiled-coil-helix domain-containing protein 4; mitochondrial intermembrane space import and assembly 40 homolog; translocase of inner mitochondrial membrane 40 hom		

For research use only.

GENE DESCRIPTION	coiled-coil-helix-coiled-coil-helix domain containing 4
GENE ID	131474
ORIGIN SPECIES	Human
ACCESSIONS	BC017082
VECTOR	pDONR223
NCBI SYMBOL	CHCHD4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctattccggcagggaagggaagatcgaatcatattgtaccaaagaagatcatgaaactccaagcagtgcagaattggtgctgatgacccaacgatccatacaggagcatggattgatactgccaatggaacattaactggaactgcccattgggggaatggcca
gctgtccctgtggagaacagtttaagtcagcctttctgcttccactatagcacggaggatcaaggggtcagactgtgtagaccagttccggccatgcaggaaatgcagaaatcccgacacctatcccaaggatgaggatgaggaggaagagagaagaagcca
gcagaacaagcagaagaacagctccattgaggccactgcaaccaagaagaggaggatcaagttaa