

RBM8A-RNA binding motif protein 8A Gene

Reference: PTXBC017088

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ribonucleoprotein RBM8A; BOV-1A; BOV-1B; BOV-1C; C1DELq21.1; DEL1q21.1; MDS014; RBM8; RBM8B; TAR; Y14; ZNRP; ZRNP1; RNA-binding protein 8A; BOV-1; RNA binding motif protein 8B; RNA-binding protein Y14; binder of OVCA1-1; ribonucleoprotein RBM8; RNA bindin		

For research use only.

GENE DESCRIPTION	RNA binding motif protein 8A
GENE ID	9939
ORIGIN SPECIES	Human
ACCESSIONS	BC017088
VECTOR	pDONR223
NCBI SYMBOL	RBM8A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggacgtgctgatcttcacgaggtctggggcggaagattcgccatggatgaggtgggacgagacattcacaaactgaaagaaaaagcgaagaaacggaagggtcgcggtttggctccgaagaggggtcccgagcgcggatgcgtgaggattatgacagcgtggagc
aggatggcgatgaacccgaccacaacgctctgttgaaggctggattctttgtaactggagtcctatgaggaagccaccgaagaagacatacacgacaattcgcagaatatggggaattaaaaacattcatctcaacctcgacagggcaacaggatattgaaggggtatacttagtt
gaatatgaaacatacaaggaagccaggctgctatggaggactcaatggccaggatttgatgggacagccatcagcgttgactgggtgtttgttgggggtccacaaaaggcaagagagaggtggccgaagacgcagcagaagtccagaccggagacgtcgctga