

TRA2A-transformer-2 alpha Gene

Reference: PTXBC017094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AWMS1; HSU53209; transformer-2 protein homolog alpha; TRA-2 alpha; TRA2-alpha; Tra2alpha; htra-2 alpha; transformer-2 alpha; transformer-2 protein homolog A; transformer 2 alpha homolog		

For research use only.

GENE DESCRIPTION	transformer-2 alpha
GENE ID	29896
ORIGIN SPECIES	Human
ACCESSIONS	BC017094
VECTOR	pDONR223
NCBI SYMBOL	TRA2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgtgtggaggaaaacttcgagggcagagagtcgtctcagtcctcaaacgggaactcctgctcgtgtaaaatcggagagcaggtcaggtctcgtatccatcaagggtttccaaactctgaatccattctcgaatcaaatccaggtcgaagtcgaag
gagacattctcatagacgttacactcgatccagatccactctcactctcatagagagcagatctcgaagtagatcatataccagaataccggcggcgaaggagccgaagccattctcaatgtctaacggagaagacatactggcagcagggcaaatccagatccaaacttgccttg
gagtgtttggcctcagttgtacacaacagagagggatctcgtgaagtattttctgatatggacattgagtggtgtcaatgtgtttatgatcagcgaactggcgatctcgaggatttctttgtgtatttgagagaatagatgactcaaggaggctatggaaagggaatggaatgga
gctggatgtagaagaattcgggtgattattctataccaagagagcgcacacccaacacccaggcatctacatggcgagaccaactcatagtggtgggggtggtggaggaggcggcggcggtggaggtggaggtggtggcagacgtcgagattcttactatgatagagatgatc
gtgggtatgacagatatgaagactatgattaccgatacagaagacgatcccttctcttattatagtcgatatagatcacgatcaagatctggtctacagccaagacgtattga