

ZNF808-zinc finger protein 808 Gene

Reference: PTXBC017109

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	zinc finger protein 808		

For research use only.

GENE DESCRIPTION	zinc finger protein 808
GENE ID	388558
ORIGIN SPECIES	Human
ACCESSIONS	BC017109
VECTOR	pDONR223
NCBI SYMBOL	ZNF808
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattgtgcaagcccttacttcacattcacacctcgtggacatcagagaatccactggacagaaatctgcaaatgcatcagtggtgcaaggtcttcagtcgaggtcactcctgcagaacatgagaaaatcattttgagataactgttccaatgcagtgagtatagcaaacatca
agcattaattgacattggagtcaattcagcattgactga