



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CREM-cAMP responsive element modulator Gene

Reference: PTXBC017117

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CREM 2beta-a protein; CREM 2alpha-b protein; CREM-2; ICER; hCREM-2; cAMP-responsive element modulator; cAMP response element modulator; inducible cAMP early repressor ICER; cAMP responsive element modulator		

For research use only.

GENE DESCRIPTION	cAMP responsive element modulator
GENE ID	1390
ORIGIN SPECIES	Human
ACCESSIONS	BC017117
VECTOR	pDONR223
NCBI SYMBOL	CREM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacacatggaacatgtagtaccacgcatgatggaagtaataacagctcttgtagacagagagcaagctgctcatgtgcagactcagactggccaaaatcaatccctgcttttagctcaggtagcagcaattgcagagacagatgaatctgcagaatcagaagggtgtaattgattctcataaac
gtagagaacatctccttcacgaagaccctcttataggaaaatactgaatgaactctcctctatgtgctctgggtgtcccaagatgaagaagagagatcagagaagaaggaagaccaccactagtattgtccatcatggcagtagccaactgacatatactgacacggggcaatacagtagtat
ctctcccaattctatagtaacagctcactcttaagctctctag

More info: Please contact us for any question: contact@proteogenix.fr