

HSPC152-hypothetical protein HSPC152 Gene

Reference: PTXBC017172

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC152; HSPC170; TRM112; TRMT11-2; hTrm112; multifunctional methyltransferase subunit TRM112-like protein; TRM112-like protein; tRNA methyltransferase 112 homolog; tRNA methyltransferase 11-2 homolog (S. cerevisiae)		

For research use only.

GENE DESCRIPTION	hypothetical protein HSPC152
GENE ID	51504
ORIGIN SPECIES	Human
ACCESSIONS	BC017172
VECTOR	pDONR223
NCBI SYMBOL	HSPC152
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactgttaccacaaatctgctgagctcgcattgtcggtgggggtgggtccctgtggttccccctgcgcctccaggccacgaggtccgtatctgccctgtggaattcaaccccaacttcgtggcgcgtatgatacctaaagtggagtggtcggcgttctggaggcggccgataacttgcgtctgatccaggtgccgaaggccggttgaggatagaggagaatgaggagtcttctgaggacatgcaccacactgctctggagggtggaagtgatagaggccacctgcagtcgccgaatctggacgtatgttccccatcagccgcggatcccaacatgctgctgagtgaaaggaaactgagagttga