



Reference: PTXBC017187

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MSK1; MSPK1; RLPK; ribosomal protein S6 kinase alpha-5; 90 kDa ribosomal protein S6 kinase 5; RSK-like protein kinase; RSKL; S6K-alpha-5; nuclear mitogen- and stress-activated protein kinase 1; ribosomal protein S6 kinase, 90kDa, polypeptide 5; ribosomal		

For research use only.

GENE DESCRIPTION	ribosomal protein S6 kinase, 90kDa, polypeptide 5
GENE ID	9252
ORIGIN SPECIES	Human
ACCESSIONS	BC017187
VECTOR	pDONR223
NCBI SYMBOL	RPS6KA5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1