

POLR2D-polymerase (RNA) II (DNA directed) polypeptide D Gene

Reference: PTXBC017205

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSRBP4; HSRPB4; RBP4; RPB16; DNA-directed RNA polymerase II subunit RPB4; DNA-directed RNA polymerase II 16 kDa polypeptide; DNA-directed RNA polymerase II subunit D; RNA polymerase II 16 kDa subunit; RNA polymerase II subunit B4; RNA polymerase II subunit		

For research use only.

GENE DESCRIPTION	polymerase (RNA) II (DNA directed) polypeptide D
GENE ID	5433
ORIGIN SPECIES	Human
ACCESSIONS	BC017205
VECTOR	pDONR223
NCBI SYMBOL	POLR2D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggtggcagcgcggtggcgcgacgtagaggaggacgctcacagctcatcttctctaaagagttgaaacagctgagacacttctaaatcagaagttcatatgctctggaacatcgaagcagcagaatgagagtcagaggacgaacaggagctctcagaagt
cttcataaaaacattaaactacagcccggttcagtcgtttcaaaaacagagaccattgccagtggttcgtagcttgcctaccagaaaaagctcataagttgagttgctgtttgccaactttgccagagactgctgaggagtcgaaggctctaatccaagcttgaggggacggtt
tgaagatgaggagctgcagcagattcttgatgataccagacaaagcgcagcttcagtattaa