



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

DKFZp761E198-DKFZp761E198 protein Gene

Reference: PTXBC017264

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

For research use only.

GENE DESCRIPTION	DKFZp761E198 protein
GENE ID	91056
ORIGIN SPECIES	Human
ACCESSIONS	BC017264
VECTOR	pDONR223
NCBI SYMBOL	DKFZp761E198
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgcagaggacaccggccctctgactgcctgagctcgggtgccatcgggtcaaggcccactccagctgttgaagctccagccggagggcctggagcccatctactctcggagctgcgtctccgtgtggaaggacagctatgacccccctggaggctgtccatgtgccctcct
gtctctggccgccctgccgccctctgtctctgcctctgcagccccgatgccggcccccgcacggctgggatgtccatgccctttacacacatccactggcttcacgtccatgccacttgcacccccctgtctgtgaacttggccagactctttctgctttcccgagctccagagggg
ggccggctgggctctggagagctcggcgattccctgcgcctccagagagctcgtgtagagctcgtgtgtgtgtcattgtccggccacaggccggcggagcttgatctccctccgccactggagctctttgtgtgtggccagcctctaccagctactgtgtagcaatccactggtcccccg

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1