



Reference: PTXBC017293

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL70; epididymis luminal protein 70; membrane-organizing extension spike protein		
For research use only.			

GENE DESCRIPTION	moesin
GENE ID	4478
ORIGIN SPECIES	Human
ACCESSIONS	BC017293
VECTOR	pDONR223
NCBI SYMBOL	MSN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]