

SESN3-sestrin 3 Gene

Reference: PTXBC017296

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SEST3; sestrin-3; sestrin 3		
For research use only.			

GENE DESCRIPTION	sestrin 3
GENE ID	143686
ORIGIN SPECIES	Human
ACCESSIONS	BC017296
VECTOR	pDONR223
NCBI SYMBOL	SESN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaccggggcgccgagcccgctgcccgcgccaactacctgctctgtaccaactgccgaaagtgtcggaaggataaaagaatcagagtgtctcaaccctgacaagaggaccaagtgcctttattccagagaaggaagttgtccaagcaaacacagtggaatgaacgtacta
actttctgtggaagaatactctacatccggtcgtctggacaacatcacacaggtcatgagtttacacactcagtaacctggagctcttcttgcggagccagttttacatgttgcgcagtgatgggtccctctctctaccatacaggcactatattgcaataatggctgcagctagacatcagtttctta
cttaataaacatgcatgtggatgaatttttaaaactggaggtattgtctgagtggttgaatggttgaatattgtccacaagactgaaaaatcttaafgaatgaataaagctgtagcacatcgacctggctgacacaaaagagcacattcagaactgtcaaaactggagaaaaaattggt
ctctgacctgaactggtacatgctgtggtcctctgacacattatcatgtttggcaagctttgttttggtagtggtatcaatccagagagagatccagaatctccaatggattcagggttaataatcagtcacaatttctgcggtttgtgatcttgctaafgacaacaacatagagaatgcattctttca
ggcagcaactttggattgtggtatcttaagttagctagagcccttaattgaaaggatgaaaagacttcaagaagaaagggaagatgaagaggcgtctcaagaagaatgagcactggtttgaaaaggagaagaagaaagtcttttgggtctctggagatactttcatttcctca
ttcaggtgctttttgcactttttgccttttag