

RPS12-ribosomal protein S12 Gene

Reference: PTXBC017321

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S12; 40S ribosomal protein S12; ribosomal protein S12		
For research use only.			

GENE DESCRIPTION	ribosomal protein S12
GENE ID	6206
ORIGIN SPECIES	Human
ACCESSIONS	BC017321
VECTOR	pDONR223
NCBI SYMBOL	RPS12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgaggaaggcattgctgctggaggtgtaattgacgttaatactgctttacaagaggttctgaagactgcctcatccacgatggcctagcacgtggaattcggaagctgccaagcccttagacaagcgccaagcccatcttgggtgcttgcacccaactgtgatgagcctatgtatg
tcaagttggtggagggcccttggctgtaacacccaataacacctaataaggttgatgacaacaagaactaggagaatgggttaggccttgtaaaattgacagagaggggaaacccgtaaaagtggttggtgagttgttagtagttaaggactatggcaaggagtcagggccaaggat
gtcattgaagagtatttcaaatgcaagaaatga