



Reference: PTXBC017338

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FUCA; tissue alpha-L-fucosidase; alpha-L-fucosidase 1; alpha-L-fucosidase I; alpha-L-fucoside fucohydrolase 1; fucosidase, alpha-L- 1, tissue		

For research use only.

GENE DESCRIPTION	fucosidase, alpha-L- 1, tissue
GENE ID	2517
ORIGIN SPECIES	Human
ACCESSIONS	BC017338
VECTOR	pDONR223
NCBI SYMBOL	FUCA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]