

RPL31-ribosomal protein L31 Gene

Reference: PTXBC017343

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L31; 60S ribosomal protein L31; ribosomal protein L31		
For research use only.			

GENE DESCRIPTION	ribosomal protein L31
GENE ID	6160
ORIGIN SPECIES	Human
ACCESSIONS	BC017343
VECTOR	pDONR223
NCBI SYMBOL	RPL31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctccgcaagaagggtggcgagaagaaaaaggccgttctgcatcaacgaagtggttaacccgagaatacaccatcaacattcacaagcgcacccatggagttgggttcaagaagcgtgcacctcgggcactcaagagattcggaaattgccatgaaggagatgggaactc
cagatgtgcgcatgacaccaggctcaacaaagctgtctggccaaaggaataaggaatgtgccataccgaatccgtgtgcggtgtccagaaaacgtaatgaggatgaagattcaccataaagctatatactttgtttacatgtacctgttaccacttcaaaaatctacagacagtcatt
gtggatgagaactaa