



Reference: PTXBC017348

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EST122234; TAPL; ATP-binding cassette sub-family B member 9; ABC transporter 9 protein; ATP-binding cassette, sub-family B (MDR/TAP), member 9; TAP-like protein; ATP binding cassette subfamily B member 9		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family B (MDR/TAP), member 9
GENE ID	23457
ORIGIN SPECIES	Human
ACCESSIONS	BC017348
VECTOR	pDONR223
NCBI SYMBOL	ABCB9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]