

ZNF765-zinc finger protein 765 Gene

Reference: PTXBC017357

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 765		
For research use only.			

GENE DESCRIPTION	zinc finger protein 765
GENE ID	91661
ORIGIN SPECIES	Human
ACCESSIONS	BC017357
VECTOR	pDONR223
NCBI SYMBOL	ZNF765
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcttcctcagggtctattgacattcaggatgtggccatagaattctctcaggaggagtggaaatgcctggaccctgctcagaggactctatacagggacgtgatgctggagaattataggaaacctggtctccctggagttgtcaggggagtgctccattggcagcacctgcctcctga