



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

HIST1H4J-histone cluster 1, H4j Gene

Reference: PTXBC017361

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H4/e; H4F2iv; H4FE; dJ160A22.2; histone H4; H4 histone family, member E; histone 1, H4j; histone cluster 1, H4j; histone cluster 1 H4 family member j		

For research use only.

GENE DESCRIPTION	histone cluster 1, H4j
GENE ID	8363
ORIGIN SPECIES	Human
ACCESSIONS	BC017361
VECTOR	pDONR223
NCBI SYMBOL	HIST1H4J
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggcgcgcgaagaagcgggaagggctgttctgcacaagggcgcgctaagcgccacgcgttaagactctgcgcacaataltcaggggcatcacaaagcgggccaatcgggccttctgtcgcgcgggcggctggaagcgcatctcgggcctcatctacgagagactcgcgggtgc
tgaagctgtctcgagaaactgtatctcggcagccctgacatctatcacagacacgcgaagcgcaagcggctaccgcgaatgtagtctgtctacgcgtctcaagcgcccaagggccgcacctctacgggttctgcgttga