

LOC492311-similar to bovine IgA regulatory protein Gene

Reference: PTXBC017422

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to bovine IgA regulatory protein
GENE ID	492311
ORIGIN SPECIES	Human
ACCESSIONS	BC017422
VECTOR	pDONR223
NCBI SYMBOL	LOC492311
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaacgcagtggtcgggctgtaataattaccataattgctgtcatgtttcccatctcagtgctgggaatcaccatgtggaaccaagcaaacgtgtgtgcatcagccggttgagttgttcaatatcaaagctga