



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MIER1-mesoderm induction early response 1 homolog (*Xenopus laevis*) Gene

Reference: PTXBC017423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MIER1 transcriptional regulator; MI-ER1; mesoderm induction early response protein 1; mesoderm induction early response 1 homolog; mesoderm induction early response 1, transcriptional regulator		

For research use only.

GENE DESCRIPTION	mesoderm induction early response 1 homolog (Xenopus laevis)
GENE ID	57708
ORIGIN SPECIES	Human
ACCESSIONS	BC017423
VECTOR	pDONR223
NCBI SYMBOL	MIER1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggaaggagaaacaaacttcagctctgaataagaagattcttcaagggaaggcgacatgccaatcatgaactctcagcccttatggttatggtatgactgttcgactacctgaagaagatgagggaaggaagaagagggaagaagaagggtgaagatgatgaagatgctgataatg
atcgaagaacagtggctgtgttggcgaataaagaaggagaataaagaagattcatcagcttcagcttggaagatgaactcctcttccaaatgacatcatgacatcactctgttgccttccagaatgccaggaataatccgccactgcgatgaaatattttgatagaatgtgaagttagaagaaga
atctgaagaagatgagaattatgctccacagaagatcgaaaaagatccgcttccactcttccagtcgaacaaatctccarttgcatttatacatactaa