

FASLG-Fas ligand (TNF superfamily, member 6) Gene

Reference: PTXBC017502

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALPS1B; APT1LG1; APTL; CD178; CD95-L; CD95L; TNFSF6; TNLG1A; tumor necrosis factor ligand superfamily member 6; CD95 ligand; Fas ligand (TNF superfamily, member 6); apoptosis (APO-1) antigen ligand 1; apoptosis antigen ligand; fas antigen ligand; mutant t		

For research use only.

GENE DESCRIPTION	Fas ligand (TNF superfamily, member 6)
GENE ID	356
ORIGIN SPECIES	Human
ACCESSIONS	BC017502
VECTOR	pDONR223
NCBI SYMBOL	FASLG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcagccctcaattaccatatacccagatctactgggtggacagcagtgccagctctccctgggcccctccagcagcagttctccctgtcccaactctgtgccagaaggcctgggtcaaggaggccaccaccaccaccgccaaccgccaactaccacctccgccgcgcgc
gccaccactgcctccactaccgctgccaccctgaagaagagagggaaccacagcacaggcctgtgtctctgtgtatgttttcatgggtctgtgcttggtaggaltgggctggggatgtttcagctcttcacctacagaaggagctggcagaactccgagagctaccagccagat
gcacacagcatcatcttggagaagcaaataggccacccagctccacccctgaaaaaaggagctgaggaagtgcccatttaacaggcaagccaactcaaggtccatgcctctggaatgggaagacacatggaattgtctgtcttctggagtgaagataagaagggtggcctt
gtgatcaatgaaactgggctgtactttgtatattcctaaagtatacttccgggtcaatcttgcaacaacctgcccctgagccacaaggtctacatgaggaacttaagatcccaggatctgggtgatgatggagggaagatgatgagctactgcactactggcagatgtggcccgccgagc
agctacctggggcagtggttcaatcttaccagtgctgatcatttatgtcaacgtatctgagctctctgtgcaatttgaggaaatcagacgtttttcggcttatataagctctaa