

TBPL1-TBP-like 1 Gene

Reference: PTXBC017559

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MGC:8389; MGC:9620; STUD; TLF; TLP; TRF2; TATA box-binding protein-like protein 1; 21-kDA TBP-like protein; TATA box-binding protein-related factor 2; TBP-like 1; TBP-related factor 2; TBP-related protein; second TBP of unique DNA protein; TATA-box bindin		

For research use only.

GENE DESCRIPTION	TBP-like 1
GENE ID	9519
ORIGIN SPECIES	Human
ACCESSIONS	BC017559
VECTOR	pDONR223
NCBI SYMBOL	TBPL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgcagacagtgatgtgcatggacattctaattacaagttagtctgtgttttagaacaagatgtcatttaactaaggagattgcttggaggagcaaatgtaattataaacgtgatgttgaaaagtattaatgaagcttagaaaacctagaattacagctacaatttggtctcag
gaaaaattattgcactggagcaacaagtgaagaagaagctaaatttggtccagacgcttagcccgtagtctgcagaaactagggtttcaggttaattacagattttaagggtgttaacgttctggcagtggtgtaacaatgccatttgaatccgtttgccagaattcacaagaacaatagacct
catgccagttacgaacctgaacttcatcctgctgtgctatcggaataatcttaagagctacattacagattttcaacaggaagatcacagtaacaggcccaatgtaaggctgttgcactgctgtggaacagatttaccatttgtgttgaaagcaggaaagaaatttataa