

NRG4-neuregulin 4 Gene

Reference: PTXBC017568

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pro-NRG4; HRG4; pro-neuregulin-4, membrane-bound isoform; neuregulin 4		
For research use only.			

GENE DESCRIPTION	neuregulin 4
GENE ID	145957
ORIGIN SPECIES	Human
ACCESSIONS	BC017568
VECTOR	pDONR223
NCBI SYMBOL	NRG4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaacagatcacgaagagccctgtgtccagtcacaagtcgtttgcctgaatgggggctttgtatgtgatactactatccagccattttaggtgcgttgaactatacaggagctcgttgaagaggttttctccaggtccagcatccaaactaaaagtaacctgttgaagcttttggcattggcggtcctagtaacacttatcatggagccttctacttcttgcaggaaaggccacttccagagagccagttcagtcagatgatatcaacctggtagagacgagcagttaccagtgtccaccacagtcattgaacaacactga