

SCRG1-scrapie responsive protein 1 Gene

Reference: PTXBC017583

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SCRG-1; scrapie-responsive protein 1; scrapie responsive gene 1; scrapie-responsive gene 1 protein; stimulator of chondrogenesis 1		

For research use only.

GENE DESCRIPTION	scrapie responsive protein 1
GENE ID	11341
ORIGIN SPECIES	Human
ACCESSIONS	BC017583
VECTOR	pDONR223
NCBI SYMBOL	SCRG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactgatgggtactgttttcaccattgggctaactttgctgctaggagttcaagccatgcctgcaaatcgctctcttgcacagaaagatactaaaagatcacaactgtcacaaccttcggaaggagtagctgacctgacacagattgatgtcaatgtccaggatcatttctaa