



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF580-zinc finger protein 580 Gene

Reference: PTXBC017698

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

zinc finger protein 580; LDL-induced EC protein

For research use only.

GENE DESCRIPTION	zinc finger protein 580
GENE ID	51157
ORIGIN SPECIES	Human
ACCESSIONS	BC017698
VECTOR	pDONR223
NCBI SYMBOL	ZNF580
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]