

C4orf18-chromosome 4 open reading frame 18 Gene

Reference: PTXBC017701

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C4orf18; AD021; AD036; protein FAM198B; expressed in nerve and epithelium during development; family with sequence similarity 198 member B		

For research use only.

GENE DESCRIPTION	chromosome 4 open reading frame 18
GENE ID	51313
ORIGIN SPECIES	Human
ACCESSIONS	BC017701
VECTOR	pDONR223
NCBI SYMBOL	C4orf18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgagggtgttgcctccacctagacagatcctgggctcaacagaccctgccgtctgtgagcaggaaagcagagttcatcaagatggcgcgccatgcccatcttttgggatgcattttcttcagcaagtaataacacccattctgttaagctcacctggggaactat
cagcagttgctgaacagaaatgctggcagaaatggccgagtagcccaagcctgaatcgggtgtactgaaatacatcatgagtggtccaagatggcactctttgatttttttacagattataatcgttagatacaaaatgctgtgattcagacctcgcaaggagatgctgtgtacaga
atggattgaggccaaaatgtgatgaccaaggttctgctgctctagcacattatccagcgaagcatgacccaaggcatttgggttttatagacaacaagggtttctttgacaggagtgagataaacttaaacctcaaatgttagaaggcatcaagagttccagcttctgcagttctgtttg
aagagccagcacttacggcagaaacttctcagctctgtttcttgataaagtgtattgggaaagtcaaggaggtagacaaggaattgaaaagcttatcatgtaataagacagagcgaataattctatcacctatatcaatgcacacgggtcaaagtattacctaataatgaatga