



Reference: PTXBC017723

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT1.4; MAGE-41; MAGE-X2; MAGE4; MAGE4A; MAGE4B; melanoma-associated antigen 4; MAGE-4 antigen; MAGE-41 antigen; MAGE-X2 antigen; cancer/testis antigen 1.4; cancer/testis antigen family 1, member 4; melanoma antigen family A, 4; melanoma antigen family A4;		

For research use only.

GENE DESCRIPTION	melanoma antigen family A, 4
GENE ID	4103
ORIGIN SPECIES	Human
ACCESSIONS	BC017723
VECTOR	pDONR223
NCBI SYMBOL	MAGEA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1