



Reference: PTXBC017734

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf163; RA410; SLY1; SLY1P; STXBP1L2; sec1 family domain-containing protein 1; syntaxin-binding protein 1-like 2; vesicle transport-related protein; sec1 family domain containing 1		

For research use only.

GENE DESCRIPTION	sec1 family domain containing 1
GENE ID	23256
ORIGIN SPECIES	Human
ACCESSIONS	BC017734
VECTOR	pDONR223
NCBI SYMBOL	SCFD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]