

PRIM2-primase, DNA, polypeptide 2 (58kDa) Gene

Reference: PTXBC017833

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRIM2A; p58; DNA primase large subunit; DNA primase 58 kDa subunit; primase, DNA, polypeptide 2 (58kDa); primase (DNA) subunit 2		

For research use only.

GENE DESCRIPTION	primase, DNA, polypeptide 2 (58kDa)
GENE ID	5558
ORIGIN SPECIES	Human
ACCESSIONS	BC017833
VECTOR	pDONR223
NCBI SYMBOL	PRIM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttttctggaagaagtggaggagctgaggttggcaggtgaccagaggagatgcttaccctcattgccttcttacttcagccacctctgaaaacatctttaataagaattgaaaacttgctattgatagagttaaattgttaaaatcagttgaaaatcttgagtgagctatg
tgaagggaactgaacaataccagagtagttagtggagagtgagcttcggaagctcaagtcttctacagagaaaaacttagaagatgaatatgaaccacgaagaagagatcatatttctatttttgcgcttgcttattgccagctgaagaacttagacgctggttcattcaacaagaatgga
tctcctcgatttagatttagtattttaccagagataaaatcaggatttcttaaggatagccaattgcagtttgaggctgtaagtataatttttag