

CH25H-cholesterol 25-hydroxylase Gene

Reference: PTXBC017843

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C25H; cholesterol 25-hydroxylase; cholesterol 25-monooxygenase; h25OH		
For research use only.			

GENE DESCRIPTION	cholesterol 25-hydroxylase
GENE ID	9023
ORIGIN SPECIES	Human
ACCESSIONS	BC017843
VECTOR	pDONR223
NCBI SYMBOL	CH25H
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctgccacaactgtccgacccccaggtcctttgagctcgggagctgttctgcagccctctgggaccacctgaggagctgggagggccctctacagtcgcccttctccgggtcatcttccatcaccacatactgggcttttgcctgccctctgtgtcctgatactgtg
ctcctgggtgcccgccttgcggcgtacaagatccacctgacttctccatccgcgcagcagctgtacttgcctggggcagacctctaccagcatgtgatgtttgtttcccgtagcgtctgtcattgggcccgcagcccgcctctgccccaggaagctccgagctgtc
ctgctgtgcaccacatcctgttctgctgtactcttcgacatggagttcttctgtgtgcacctgtgcaccacaagggtccctggctgtacgcacacctccacaagggtgcaccaccagaactgtcctctgttcgctgggaacgagtatatgagcgtctgggaactgtttcttgggcttct
tcgacatgatgaacgtcacactgtcgggtgcccccgtcaccacctgacctccacgtgggtcaacatctggttccgtgggagaccactccggtacaaacttccctgggtccactcacagactgtgccctcgggtgttacgggggtgtgtgcaccagacctgcatcacttcaact
ttaactgaacttgcctccgtactttacacactgggacaaaatactgggaacgtgcggactgcatctgtccagcgcgtga