

TFF3-trefoil factor 3 (intestinal) Gene

Reference: PTXBC017859

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ITF; P1B; TFI; trefoil factor 3; polypeptide P1.B; trefoil factor 3 (intestinal)		
For research use only.			

GENE DESCRIPTION	trefoil factor 3 (intestinal)
GENE ID	7033
ORIGIN SPECIES	Human
ACCESSIONS	BC017859
VECTOR	pDONR223
NCBI SYMBOL	TFF3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctggggctggtctctgctgctccagctctgctgaggagtacgtggcctgtctgcaaacagtggtgcccgtgccagccaaggacagggtggactgcggctacccccatgtcaccccaaggagtgcacaacacggggctgctgctttgactccaggatccctggagtgct
tgggtttcaagcccctgcaggaagcagaatgcaccttctga