

FBXO36-F-box protein 36 Gene

Reference: PTXBC017869

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Fbx36; F-box only protein 36; F-box protein 36		

For research use only.

GENE DESCRIPTION	F-box protein 36
GENE ID	130888
ORIGIN SPECIES	Human
ACCESSIONS	BC017869
VECTOR	pDONR223
NCBI SYMBOL	FBXO36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtcgtggcgtccggagactctttgaaactgtaggacaaggccgcctagcaaaagactattaccagttactgtcacccggtctcaggtaatctttagatgggtgaagatctcttaaggagtgtatcgcatacaaaaacctggagaagcaaaagaacccatgaagacttcctagagaattcacatctcaaggttaaggacgaacatattatataccctgttaagtctcattagttagtgtagattttgttaa