

P2RY12-purinergic receptor P2Y, G-protein coupled, 12 Gene

Reference: PTXBC017898

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	purinergic receptor P2RY12; ADPG-R; BDPLT8; HORK3; P2T(AC); P2Y(12)R; P2Y(AC); P2Y(ADP); P2Y(cyc); SP1999; P2Y purinoceptor 12; ADP-glucose receptor; G-protein coupled receptor SP1999; Gi-coupled ADP receptor HORK3; P2Y12 platelet ADP receptor; purinergic		

For research use only.

GENE DESCRIPTION	purinergic receptor P2Y, G-protein coupled, 12
GENE ID	64805
ORIGIN SPECIES	Human
ACCESSIONS	BC017898
VECTOR	pDONR223
NCBI SYMBOL	P2RY12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaagcgtcgacaatcaccctctgcgcctgggaacaccagctctgtgaccagagactacaaaatcaccaggctctctccactgctctacactgtcctgtttttgttgacttatcacaatggcctggcgatgaggattttcttcaaatccggagtaaatcaaatcttatattttcttaa
gaacacagtcatttctgatcttctcatgattctgactttccattcaaaattcttagtgatgccaactgggaacaggaccactgagaactttgtgtgtcaagttaacctccgctcatattttattcacaatgtatcatgattttcattctgggactgataactatcgaatcctaccagaagaccaccagg
ccattttaaactccaaccccaaaatctcttgggggctaagattctctctgtgtgctatctgggcaattcatgttcttactctcttgcctaactgattctgaccaacaggcagccgagagacaagaatgtgaagaatgctcttcttcaatcagatctggctagctggcatgaaatagtaaat
acatctgcaagtcattttctgattaatttctaattgtattgtatgtataactcattacaaaagaactgtaccggctacagtaagaacgaggggtgtaggtaaagtcgccaggaagggtaacgtcaaaagtgttcaattatcattgctgtattcttattgtttgtcttcttcttccatttgcggaaat
ccttacacctgagccaaacccgggatgtcttctgactgactgctgaaaatactctgttctatgtgaaagagagcactctgtgtgtaactcttaaatgcatgcctggatccgttcatctatttttcttgcagtccttcagaaatctctgataagatgtggaagtgcctcaattctgcaacatct
ctgtccaggacaataggaaaaagaacaggatgggtgtgacccaaatgaagagactccaatgtaa