

CPLX3-complexin 3 Gene

Reference: PTXBC018026

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPX-III; CPXIII; Nbla11589; complexin-3; CPX III; complexin III; complexin 3		
For research use only.			

GENE DESCRIPTION	complexin 3
GENE ID	594855
ORIGIN SPECIES	Human
ACCESSIONS	BC018026
VECTOR	pDONR223
NCBI SYMBOL	CPLX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgttcattggtgaagaccatggtggcgccagctgaagaacctcactgggagcctgggagcgccgaggataaggagatgggacaagtggcagccgaagctcagggcatgagccggaggagtagcaggatgatcagaagcaactcgtggaagagaagatggagc
gggatgcacagttcacacagaggaaggcagagcgccacactgcggagccacttcgagacaataccggctacccaagaacgagacagatgagagccagatccagatggcaggtggagacgtggagctgccccgggagctggccaagatgacgaggagacacagagga
ggaggaggagaaggcctcagtcctggcgagctggccagccttcctggctgaactgggctcactcaaggacaaggccagggccacactggggatctcaagcaatcagctgagaagtgtcacgtcatgtga