

## WIF1-WNT inhibitory factor 1 Gene

Reference: PTXBC018037

|                               |                                |                     |                |
|-------------------------------|--------------------------------|---------------------|----------------|
| Product type:                 | DNA                            | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | WIF-1; wnt inhibitory factor 1 |                     |                |
| <b>For research use only.</b> |                                |                     |                |

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | WNT inhibitory factor 1                                                                           |
| <b>GENE ID</b>          | 11197                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC018037                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | WIF1                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcccggaggagcgcttcctgccgcgcgtctggtctggagcatcctctgtgctgctggcactcgggcgaggcgccggcgccgagggagagcctgtacatgatgatgctcaccaggcaagagtagctataggattgaagaagatcctgattgttcgga  
ggggaaatggcacttttacacatgatttcagaaaagcgaacagagaaatgccagctattctgtcaatatccattccatgaattttacctggcaagctgcaggcgaggcagaatactctatgaattcctgctctgctccctggataaaggcatcatggcagatccaacgtcaatgtccc  
tctgctgggaacagtgctcacaaggcatcagttgtcaagttggttcccatgtcttgaaaaacaggatggggtggcagcattgaagtgatgtgattgtatgaattctgaaggcaacaccattctcaaacacctcaaatgctatctctttaaacaatgtcaacaagctgagtgccaggc  
gggtgccgaaatggaggctttgtaatgaaagacgcatcgcgagtgctgctggttccacggacctcactgtgagaagccctttgtacccacgagtgatgaatgggtgactttgtgactcctggtttctgcatctgcccacctggattctatggagtgaactgtgacaaagcaaatg  
ctcaaccactgctttaatggaggacgtgtttcaccctggaaaatgtatttgccctccaggactagaggagagcagtggtgaaatcagcaaatgccacaacctgtcgaatggaggtaaatgcattgttaaaagcaaatgtaagtgtccaaaggtaccaggagacctctgttcaaa  
gcctgtctgcgagcctggtgtgtgcacatggaacctgccatgaacccaacaatgccaatgtcaagaaggttgcatggaagacactgcaataaagggtacgaagccagcctacatgacctgaggccagcaggcgcccagctcaggcagcacacgcttcaaaaaggcc  
gaggagcgccgggacccctgaatccaattacatctgtga