

TAC1-tachykinin, precursor 1 Gene

Reference: PTXBC018047

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Hs.2563; NK2; NKNA; NPK; TAC2; protachykinin-1; PPT; neurokinin 1; neurokinin 2; neurokinin A; neurokinin alpha; neuromedin L; neuropeptide K; neuropeptide gamma; preprotachykinin; protachykinin; substance K; tachykinin 2; tachykinin, precursor 1 (substan		

For research use only.

GENE DESCRIPTION	tachykinin, precursor 1
GENE ID	6863
ORIGIN SPECIES	Human
ACCESSIONS	BC018047
VECTOR	pDONR223
NCBI SYMBOL	TAC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaatcctcgtggccttggcagctcttttctgtctccactcagctgtttgcagaagaataggagccaatgatgatctgaattactgtccgactggtagcagcgcaccagatcaaggagggaactgccggagccctttgagcatcttctgcagagaatcgcccgagaccaagcctc
agcagttcttggattaatgggcaaacgggatgctgattctcaattgaaaaacaagtggccctgttaaggctcttattggacatggccagatctctcacaagaacataaaacagattcttgggactaatgggcaaaagagctttaattctgtgcttatgaaaggagtgcaatgcaga
attatgaaagaagacgttaa