

GAGE2E-G antigen 2E Gene

Reference: PTXBC018052

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GAGE-2E; GAGE8; G antigen 2E		

For research use only.

GENE DESCRIPTION	G antigen 2E
GENE ID	26749
ORIGIN SPECIES	Human
ACCESSIONS	BC018052
VECTOR	pDONR223
NCBI SYMBOL	GAGE2E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggcgagggaagatcgacatcggcctagaccaagacgctacgtagagcctcctgaaatgattggcctatcgcccgagcagttcagtgatgaagtgggaaccagcaacacctgaagaaggggaaccagcaactcaacgtcaggatcctgcagctgctcaggaggagag
gatgaggagacatctgcaggtcaaggccgaagcctgaagctgatagccaggaacaggggtcacccacagactgggtgtgagtgtaagatggtcctgatggcgaggatggaccgccaaatccagagggtgaaaacgcctgaagaaggtgaaaagcaatcacagtgttaa