

RIT2-Ras-like without CAAX 2 Gene

Reference: PTXBC018060

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GTP-binding protein Rit2; RIBA; RIN; ROC2; GTP-binding protein Roc2; Ric-like, expressed in neurons; ras-like protein expressed in neurons; ras-like without CAAX protein 2; Ras like without CAAX 2		

For research use only.

GENE DESCRIPTION	Ras-like without CAAX 2
GENE ID	6014
ORIGIN SPECIES	Human
ACCESSIONS	BC018060
VECTOR	pDONR223
NCBI SYMBOL	RIT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggtagaaatgaagccagctgctccccgggcagcgcacagcgggtccagagagtagaaggtggtatgctgggagcagggggagtggtaaaagcgaatgacaatgcagtttattagtcacagttccctgattatcatgacccactatagaagatgcttataagaccagg
tcaggattgacaatgagccagcttacttggacatcttggacactgtggccaggcagaattcacagccatgcgggagcagtagatgcaggtggggaaggcttcatctctactccgtcactgacgtcaatcatttcaggaggctgccaagtttaagagctcattttcaggtccgcca
cacctatgaaatccccctgtgtgtgtggtgtaacaaattgatctggaacagttccgccaggtttctacagaagaaggttgagttcttgcgaagaataatgtgtgtttttttgagacctctgcagccctcagattctgtattgatgatgcttttcatggcttagtggggaattcgcaagaagga
gtccatgccatccttgatggaaaagaactgaagagaaaagactgcctgtggaagaagctcaaaggttcttgaagaagaagagagaaaatatgacatga