

RGS3-regulator of G-protein signaling 3 Gene

Reference: PTXBC018072

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2PA; RGP3; regulator of G-protein signaling 3; regulator of G-protein signalling 3		
For research use only.			

GENE DESCRIPTION	regulator of G-protein signaling 3
GENE ID	5998
ORIGIN SPECIES	Human
ACCESSIONS	BC018072
VECTOR	pDONR223
NCBI SYMBOL	RGS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgagggcatgtacctactcgcaacgggaacctgcagagggcgacacacgatgaagggaagccaaggacatgaagaacaagctggggatcttcagacggcggagtgagtcacctggagccctcccggggcaaggcagacaaaatgatgaagtcattcaagcccactcag
agggaagccctcaagtggggcgagtccttggaagctgctgttcacaatacgggttagcagtggttccaagccttcttcgactgagttcagtgaggagaaatctggagttctgttggtgctgttgaggacttcaagagggtcaagtcacagtcacagatggcatccaaggccaagaagatc
tttgcataatcgcgcgacccagcagcgaagggtcaacctggactcctacacggggagcacaccaaggacaacctgcagagcgtcacgcgggctgcttcgacctggcacagaagcgcacatctcgggctcatgaaaaggactcgtacctcgtttctccgttctgacctctac
ctggaccttattaaccagaagaagatgagtcctcccgctttag