



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM170A-transmembrane protein 170A Gene

Reference: PTXBC018082

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM170; transmembrane protein 170A; transmembrane protein 170		
For research use only.			

GENE DESCRIPTION	transmembrane protein 170A
GENE ID	124491
ORIGIN SPECIES	Human
ACCESSIONS	BC018082
VECTOR	pDONR223
NCBI SYMBOL	TMEM170A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcgcgcaggggagcggcggcagcggcgggtgcggcgggtctctgcagcagatctcgcagctgaagttgtgccgcgggtgggcaacgggacctgtgcccaactctacttccctctgctcttccagagatgtggtatggtgtattctctgtgggacctggtgtctctctcttcttctatctcctctgctgtaattcgtgcccttaccctcagacacacaaatattgtaggtctatctgtctgaacatctctgagcggcatctggggaccaattactctggaaactctcaagtgctcagctattgctggaggtttaccgagcagcaggggaaggaaatgatacattggaagccctcacatgggcactgtgacagacattttgcgtctgtgtgcctctcttttaccgatttagctactctatag