

RPS21-ribosomal protein S21 Gene

Reference: PTXBC018140

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HLDF; S21; 40S ribosomal protein S21; 8.2 kDa differentiation factor; human leukemia differentiation factor; ribosomal protein S21		

For research use only.

GENE DESCRIPTION	ribosomal protein S21
GENE ID	6227
ORIGIN SPECIES	Human
ACCESSIONS	BC018140
VECTOR	pDONR223
NCBI SYMBOL	RPS21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagacgacgccgagtcgtggacctgtacgtgccgcggaatgctccgtagcaatcgcacatcgggtccaaggaccacgcacatccagatgaacgtggccgaggttgacaaggtcacaggcaggtttaatggccagtttaaaacttatgctatctcggggccattcgtaggatgggtgagtcagatgattcattctccgattggccaaggccgatggcatcgtctcaagtaa